



Computer Specification for Pyris Software

3.0 GHz Intel® Core 2 Duo, 1333 MHz Front Side Bus (FSB), 6M L2 cache

Chassis: Tower

Factory Installed Operating System: Microsoft® Windows® XP Professional SP3

Memory: 2 GB, Non-ECC, 1066 MHz, DDR3, 2 x 1GB, Four DIMM slots

Hard Drive: 160 GB SATA, 7200 RPM

Networking: Integrated Gigabit Ethernet

External I/O Ports: 8 USB 2.0 (2 front, 6 rear), 1 Ethernet (RJ45), 2 serial (9-pin), 1 VGA (DB-15, Display Port) out

Expansion Slots:

- Slot 1: half-length, full-height, PCIe 2.0 x16 (75w max)
- Slot 2: half-length, full-height, PCIe x1
- Slot 3: half-length, full-height, 32-bit PCI 2.3
- Slot 4: half-length, full-height, 32-bit PCI 2.3

Removable Media: DVD+/-RW SATA

Video: Integrated video, Intel® GMA4500

Floppy Drive 3.5" 1.44MB diskette drive standard

Mouse USB optical mouse with scroll

Keyboard: USB Keyboard

Audio: Integrated High Definition Audio